

GLOSTAVENT® HELIX
SPARE PARTS LIST

> Glostavent® Helix spare parts

CATEGORY	PART NAME	PAGE	CATEGORY	PART NAME	PAGE
Autoclavable	Silicone adult autoclavable tubing	4	Concentrator	Equalization valve 220v - 8lpm	14
Circuits	Silicone paediatric autoclavable tubing	4		Muffler - 10lpm	14
	Patient y piece	4		Muffler - 8lpm	14
Oxygen	Straight Connector 22F - 15M	5		Sieve beds (pair) - 10lpm	14
Circuits	Oxygen mask, adult	5		Sieve beds (pair) - 8lpm	14
	Flextube 15mm	5		Compressor assy 220v 50hz - 10lpm	15
Patient	0.5 Litre reservoir bag	6		Compressor assy 220v 50hz - 8lpm	15
Circuit	1 litre reservoir bag	6		Capacitor - 10lpm	15
	2 Litre reservoir bag	6		Capacitor - 8lpm	15
	Scavenging tubing 30mm	6	Vaporiser	Isoflurane vaporiser, 150 ml	16
	Straight connector 22f - 15m	6		Sevoflurane vaporiser, 150 ml	16
	Ayres t piece paediatric circuit	6		Mini funnel for vaporiser	16
	Adult mask	7		Vaporiser service kit	16
	Paediatric mask	7	UPS	Set of replacement batteries - B400 UPS	17
	Straight connector 22m x 22m	7		Set of replacement batteries - C400 UPS	17
	Tee connector 22m x 22m x 22f	7		Replacement ups mains lead UK plug	17
	Straight connector 30m x 30m	7		Replacement ups mains lead EU plug	17
	Straight connector 30m x 30f	7		Replacement certaups	17
	One way valve 22m x 22f	8	General	8003 series WOB-L piston pump	18
	Reusable y-piece 22m-22m-22m/15f	8		12 Volt push button	18
	Paediatric self-inflating bag	8		12 Volt toggle switch	18
	Adult self-inflating bag	8		Water trap	18
	Double swivel elbow 15m-22m/15f	8		Battery charger	18
	Silicone adult autoclavable tubing	9		6mm clear length tubing 1m	18
	Silicone paediatric autoclavable tubing	9		6mm black length tubing 1m	19
Helix	Pressure gauge	10		Over pressure valve	19
Ventilator	Circuit board for ventilator	10		Inspiratory valve	19
	Loom for ventilator circuit board	10		Expiratory valve	19
	Solenoid valve for gv ventilator	10		Ventilator flow regulator	19
	Black finish aluminium knob, 22mm	10		Air regulator	20
	Black finish aluminium knob, 28mm	10		Air flowmeter	20
	Fuse for ventilator circuit board	11		Oxygen flowmeter	20
	Control panel battery 12v 10 Ah	11		Scavenger unit	20
Concentrator	10lpm oxygen concentrator	12		Flowmeter flow regulator	20
	Final bacteria filter	12		Non return valve, 6mm	20
	Foam filter	12		Non return valve, 4mm	21
	Filter felt intake	12		Oxygen flush valve	21
	Compressor service kit	12		Elbow block 34mm	21
	Heat exchanger modified	12		Offset elbow block 34mm	21
	Circuit board for 10lpm	13		Bullnose/Pin Index Adaptor	21
	Circuit board upgrade kit	13		Oxygen cylinder bullnose pressure regulator	21
	O2 monitor board for or 10lpm	13		Pressure regulator	22
	New life intensity single flow Concentrator (3 led)	13		Metal tray	22
	O2 monitor board for 8lpm and 10lpm	13		Oxygen hose 2m	22
	Newlife intensity concentrator (1 led)	13		Male connector	22
	Solenoid valve 220v	13		Female connector	22
	Valve stem rebuild kit	13		22M x 22F black block connector	22
	Equalization valve 220v - 10lpm	14			

CATEGORY	PART NAME	PAGE
General	Disposable peep valve	23
	Peep valve casing	23
	Reusable AMBU peep valve 20	23
	Peep valve casing	23
	Senseca o2 analyser/ oxygen measuring device	23
	Greisinger GOEL 370 Oxygen Sensor	23
	Fuse of ventilator circuit board	24
	AVS 220v 50 hz plus cable	24
	O2 probe	24
	Check valve	24
	Circuit holder	24

Contact information

Diamedica (UK) Ltd

Unit 1 Kingdom Fields,
Bratton Fleming,
Devon, UK
EX31 4EN

Tel: +44 (0)1598 710066

Email: info@diamedica.co.uk


Web: www.diamedica.co.uk

Technical Support:

Email: support@diamedica.co.uk

WhatsApp: +44 (0)7716 503156

> Autoclavable circuits spare parts



Reference Code	002-001-1150
Description	Silicone adult autoclavable tubing
Quantity	2 / per patient per set-up

Reference Code	002-001-2150
Description	Silicone paediatric autoclavable tubing
Quantity	2 / per patient per set-up



Reference Code	002-011-1002
Description	Patient Y piece
Quantity	1 / per patient per set-up


Oxygen circuits spare parts <



Reference Code	1962000
Description	Straight Connector 22F - 15m
Quantity	2 / per patient per set-up



Reference Code	1188015
Description	Oxygen mask, adult
Quantity	1 / per patient per set-up



Reference Code	1576000
Description	Flextube 15mm
Quantity	1 / per patient per set-up

> Patient circuit spare parts




Reference Code	2805000
Description	0.5 litre reservoir bag
Quantity	1



Reference Code	2810000
Description	1 litre reservoir bag
Quantity	1




Reference Code	2820000
Description	2 litre reservoir bag
Quantity	1



Reference Code	1570 000 B
Description	Scavenging tubing 30mm, 3m length
Quantity	1



Reference Code	1962000
Description	Straight connector 22F - 15M
Quantity	1



Reference Code	2120000
Description	Ayres T piece paediatric circuit
Quantity	1

Patient circuit spare parts <



Reference Code	1516111
Description	Adult mask
Quantity	1



Reference Code	1514111
Description	Paediatric mask
Quantity	1



Reference Code	1960000
Description	Straight connector 22M x 22M
Quantity	1



Reference Code	1982000
Description	Tee connector 22M x 22M x 22F
Quantity	1



Reference Code	1970001
Description	Straight connector 30M x 30M
Quantity	1



Reference Code	1971001
Description	Straight connector 30M x 30F
Quantity	1

> Patient circuit spare parts



Reference Code	1921000
Description	One way valve 22M x 22F
Quantity	1



Reference Code	002-011-1002
Description	Reusable Y-piece 22M-22M-22M/15F
Quantity	1



Reference Code	038-72-840-H
Description	Paediatric self-inflating bag
Quantity	1



Reference Code	038-71-840-H
Description	Adult self-inflating bag
Quantity	1



Reference Code	1897000
Description	Double swivel elbow 15M-22M/15F
Quantity	1



Reference Code	002-001-1150
Description	Silicone adult autoclavable tubing
Quantity	2 / per patient per set-up



Reference Code	002-001-2150
Description	Silicone paediatric autoclavable tubing
Quantity	2 / per patient per set-up


> Helix ventilator spare parts




Reference Code	86328427
Description	Pressure gauge
Quantity	1



Reference Code	CBV-GV
Description	Circuit board for ventilator
Quantity	1



Reference Code	CBV-L
Description	Loom for ventilator circuit board
Quantity	1



Reference Code	SY3220-6LOU-M5-Q
Description	Solenoid valve GV Helix
Quantity	1



Reference Code	498-839
Description	Black finish aluminium knob, 22mm
Quantity	1



Reference Code	498-823
Description	Black finish aluminium knob, 28mm
Quantity	1



Reference Code	563-609 A
Description	Fuse for ventilator circuit board
Quantity	1

Reference Code	660-3627
Description	Control panel battery 12v 10 Ah
Quantity	1


> Concentrator spare parts



Reference Code	AS099-10-MGV
Description	10lpm oxygen concentrator Configured for use with Glostavent (also replacement for units with 8lpm concentrators)
Quantity	1



Reference Code	PV5LD-651
Description	Final bacteria filter * For both 8lpm and 10lpm Newlife Intensity oxygen concentrator
Quantity	1



Reference Code	FI002-1
Description	Foam filter *For both 8lpm and 10lpm Newlife Intensity oxygen concentrator
Quantity	1



Reference Code	MI161-2
Description	Filter felt intake For 10lpm Newlife Intensity Oxygen Concentrator *Only applicable to 10lpm concentrator
Quantity	1



Reference Code	CO012-11
Description	Compressor service kit *For 8lpm and 10lpm oxygen concentrator
Quantity	1




Reference Code	HX002-2-MOD
Description	Heat exchanger modified *For 8lpm and 10lpm Oxygen concentrator
Quantity	1

Concentrator spare parts <



Reference Code	CB200-2S
Description	Circuit board for 10lpm 220v New life Intensity single flow concentrator (3 LED) *Only for 10lpm concentrators with 3 LED **For 10lpm concentrators with 1 LED refer to KI608-2**
Quantity	1



Reference Code	KI608-2
Description	Circuit board upgrade kit CB200-2 NewLife 10L (230 V, 50 Hz). *Only for 10lpm concentrator with 1 LED
Quantity	1



Reference Code	CB196-1
Description	O2 monitor board for or 10lpm New life Intensity single flow concentrator (3 LED) *Only for 10lpm concentrators with 3 LED **For 10lpm concentrators with 1 LED refer to CB100-1**
Quantity	1



Reference Code	CB100-1
Description	O2 monitor board for 8lpm and 10lpm Newlife Intensity concentrator (1 LED) * Only for 8lpm concentrators and 10lpm concentrators with 1 LED
Quantity	1



Reference Code	VA001-2
Description	Solenoid valve 220v *For 8lpm and 10lpm Newlife Intensity oxygen concentrator
Quantity	1



Reference Code	VA034-1
Description	Valve stem rebuild kit *For 8lpm and 10lpm Newlife Intensity oxygen concentrator
Quantity	1

> Concentrator spare parts



Reference Code	VA003-2
Description	Equalization valve 220v For Newlife Intensity 10lpm concentrator
Quantity	1



Reference Code	VA003-10
Description	Equalization valve 220v For Newlife Newlife Intensity 8lpm Concentrator
Quantity	1



Reference Code	MU079-1S
Description	Muffler For 10lpm Newlife Intensity Concentrator
Quantity	1



Reference Code	MU0432S
Description	Muffler For 8lpm Newlife Intensity Concentrator
Quantity	1



Reference Code	BE186-2S
Description	Sieve beds (pair) For 10lpm Newlife Intensity oxygen concentrator
Quantity	1



Reference Code	BE157-1S
Description	Sieve beds (pair) For 8lpm Newlife Intensity oxygen concentrator
Quantity	1

Oxygen Concentrator spare parts <



Reference Code	CO306-2S
Description	Compressor assy 220v 50Hz For Newlife Intensity 10lpm concentrator
Quantity	1



Reference Code	CO006-19S
Description	Compressor assy 220v 50Hz For Newlife Intensity 8lpm concentrator
Quantity	1



Reference Code	CC012-9
Description	Capacitor For Newlife Intensity 8lpm Concentrator
Quantity	1



Reference Code	CC012-7
Description	Capacitor For Newlife Intensity 10lpm concentrator
Quantity	1

> Vaporiser spare parts




Reference Code	VAP/I
Description	Isoflurane vaporiser, 150 ml
Quantity	1



Reference Code	VAP/S
Description	Sevoflurane vaporiser, 150 ml
Quantity	1



Reference Code	5060021831255
Description	Mini funnel for vaporiser
Quantity	1



Reference Code	3400
Description	Vaporiser service kit
Quantity	1



Reference Code	RBS/UPS (Compatible with models with B400)
Description	Set of replacement batteries For B400 UPS
Quantity	1




Reference Code	RBS/UPS-C400 (Compatible with models with C400)
Description	Set of replacement batteries For C400 UPS
Quantity	1



Reference Code	42 0542
Description	Replacement UPS mains lead UK plug
Quantity	1



Reference Code	42 0563
Description	Replacement UPS mains lead EU plug
Quantity	1




Reference Code	C400-020-B12
Description	Replacement CertaUPS
Quantity	1

> General spare parts



Reference Code	80030055
Description	8003 Series WOB-L piston pump 12v DC
Quantity	1



Reference Code	1001
Description	12 volt push button
Quantity	1



Reference Code	1002
Description	12 volt toggle switch
Quantity	1




Reference Code	398 9091-MOD
Description	Water trap
Quantity	1



Reference Code	1204138
Description	Battery charger
Quantity	1



Reference Code	686-2680-A
Description	6mm clear length tubing 1m
Quantity	1



Reference Code	721-4007-A
Description	6mm black length tubing 1m
Quantity	1



Reference Code	OPV/001
Description	Over pressure valve
Quantity	1



Reference Code	PCV IN
Description	Inspiratory valve
Quantity	1



Reference Code	PCV EX
Description	Expiratory valve
Quantity	1




Reference Code	748-0709
Description	Ventilator flow regulator
Quantity	1



Reference Code	AR20-F01B-B
Description	External air regulator
Quantity	1

> General spare parts



Reference Code	435-9609-MOD
Description	Air regulator
Quantity	1



Reference Code	FL-2040-A
Description	Air flowmeter
Quantity	1



Reference Code	FL-2040-O
Description	Oxygen flowmeter
Quantity	1



Reference Code	SCAV-GV
Description	Scavenger unit
Quantity	1



Reference Code	492-8570
Description	Flowmeter flow regulator
Quantity	1



Reference Code	367-0624
Description	Non return valve, 6mm
Quantity	1



Reference Code	367-0618
Description	Non return valve, 4mm
Quantity	1



Reference Code	885-2095
Description	Oxygen flush valve
Quantity	1



Reference Code	CP002
Description	Elbow block 34mm
Quantity	1



Reference Code	CP002W
Description	Offset elbow block 34mm
Quantity	1



Reference Code	160201-02
Description	Bullnose/Pin Index Adaptor
Quantity	1




Reference Code	819-0017
Description	Oxygen cylinder bullnose pressure regulator
Quantity	1

> General spare parts



Reference Code	50.780.00.10.00
Description	Pressure regulator
Quantity	1



Reference Code	GN14-65
Description	Metal tray
Quantity	1




Reference Code	OXY-A/SH-A
Description	Oxygen hose 2m
Quantity	1



Reference Code	M Connector
Description	Male connector
Quantity	1



Reference Code	F Connector
Description	Female connector
Quantity	1



Reference Code	CP003
Description	22M x 22F black block connector
Quantity	1



Reference Code	199103001
Description	Disposable PEEP valve
Quantity	1



Reference Code	199-PEEP-CAS
Description	PEEP valve casing For 199103001 PEEP Valve
Quantity	1



Autoclavable

Reference Code	A000213000
Description	Reusable Ambu peep valve 20
Quantity	1



Reference Code	A213-PEEP-CAS
Description	PEEP valve casing For A000213000 PEEP Valve
Quantity	1



Reference Code	ECO410-35-CAN
Description	Senseca O2 analyser/ oxygen measuring device
Quantity	1



Reference Code	60-3737
Description	Greisinger GOEL 370 Oxygen Sensor
Quantity	1

> General spare parts



Reference Code	563-609 A
Description	Fuse of ventilator circuit board
Quantity	1



Reference Code	91300000-MOD
Description	AVS 220v 50 hz plus cable
Quantity	1



Reference Code	RP1
Description	O2 probe
Quantity	1



Reference Code	CH- VALVE001
Description	Check valve
Quantity	1



Reference Code	CH001
Description	Circuit holder
Quantity	1