

Oxygen Flow Splitter

Specifications

Product Name	Flow Splitter
Flow Meter Calibration	0-2LPM/0-5LPM adjustable
Accuracy	0.125LPM
Size	270*138*160mm
Inlet Pressure	48-345Kpa
Operating Temperature	5-45°C
Storage Temperature	1-60°C
Accuracy	±5°C
Usage	Connect with oxygen concentrator
OEM&ODM	Acceptable
Accessories	1pc connection tubing-1.5m 5pcs Cannula(neonatal)-2m each

