

BABY CPAP 10
SPARE PARTS LIST

> Baby CPAP 20 spare parts

CATEGORY	PART NAME	PAGE	CATEGORY	PART NAME	PAGE
Autoclavable circuits	Silicone autoclavable tubing	3	Water Trap (external) circuits	Elastomeric connector	10
	Straight connector	3		Straight connector	10
	Elastomeric connector	3		Fixed elbow - 15M-22M/15F	10
NON - Autoclavable circuits	Flextube	4		Water trap	10
	Elastomeric connector	4	Oxygen Concentrator circuits	Regulator	11
	Connector	4		Acrylic Flow Meter	11
	CPAP pipe fitting white acetal	4		Door intake filter package	11
Patient circuit	Straight connector	4		Control circuit board	11
	Humidifier bottle	5		Control circuit board starting with serial number B	11
	Bubble PAP valve bottle	5		Sieve beds	11
	Baby bonnet XS - size 1	5		T25 screwdriver	12
	Baby bonnet S - size 2	5		Water trap (Integral)	12
	Baby bonnet M - size 3	5		Straight connector - 6mm	12
	Baby bonnet L - size 4	5		Cabinet Air Filter	12
	Baby bonnet XL - size 5	6		Black length tubing	12
	Silicone nasal prong - size 0	6		Clear length tubing	12
	Silicone nasal prong - size 1	6		Flow Controller	13
	Silicone nasal prong - size 2	6		Push in straight fitting	13
	Silicone nasal prong - size 3	6		Voltage protector	13
	Silicone nasal prong - size 4	6		Exhaust muffler	13
	Silicone nasal prong - size 5	7		UK plug (Connector Converter)	13
	Silicone nasal prong - size 6	7		EU (Power Lead)	13
	Silicone nasal prong - size 7	7		UK (Power Lead)	14
	Nasal Cannula kit - Small	7		Power switch	14
	Nasal Cannula kit - Medium	7		Intake Bacteria Filter	14
	RAM Cannula Micro Preemie - size 0	8		Final Bacteria Filter	14
	RAM Cannula Preemie - size 1	8		Circuit Breaker, 3A	14
	RAM Cannula Newborn - size 2	8		Castor, non-locking	14
	RAM Cannula Infant - size 3	8		Castor, Locking	15
	RAM Cannula Small Paediatric - size 4	8		Capacitor, 220V AC	15
	RAM Cannula Medium Paediatric - size 5	8		European compressor	15
	RAM Cannula Large Paediatric	9		Compressor Filter	15
	Y-piece for RAM Cannula	9		Motor Mounts	15
	Straight connector 15F - 10M	9		Compressor Rebuild Set	15
				Cooling Fan	16
				Intake Canister	16
				Sieve Beds Check Valves	16
				4-Way Rotary Valve	16
				Accumulator Tank	16
				Pressure Regulator	16

Contact information

Diamedica (UK) Ltd

Unit 1 Kingdom Fields,
Bratton Fleming,
Devon, UK
EX31 4EN

Tel: +44 (0)1598 710066

Email: info@diamedica.co.uk

Web: www.diamedica.co.uk

Technical Support:

Email: support@diamedica.co.uk

WhatsApp: +44 (0)7716 503156

Option one – Autoclavable Patient Limbs <



Reference Code	002-002-3150E
Description	Silicone autoclavable tubing - 10F cuff, length 150cm, reusable
Quantity	2 / per patient per set-up

Reference Code	1661000
Description	Straight connector - 15M/11.5 - 10M
Quantity	2 / per patient per set-up



Reference Code	1702000
Description	Elastomeric connector - 15F - 6-9mm spigot
Quantity	1 / per patient per set-up

> Option two – Single-use Patient Limbs



Reference Code	1576000
Description	Flextube - 15mm
Quantity	2 / per patient per set-up



Reference Code	1702000
Description	Elastomeric connector 15F - 6-9mm spigot
Quantity	1 / per patient per set-up



Reference Code	1945000
Description	Connector - 15m - 15m
Quantity	1 / per patient per set-up



Reference Code	PF-CPAP
Description	CPAP pipe fitting white acetal
Quantity	2 / per patient per set-up



Reference Code	1962000
Description	Straight connector - 22F - 15M
Quantity	1 / per patient per set-up



Reference Code	02-140
Description	7600-0 Humidifier bottle
Quantity	1

Reference Code	821135
Description	Reusable humidifier bottle
Quantity	1



Reference Code	BC01/AN0018
Description	Bubble PAP valve bottle
Quantity	1

Reference Code	AN0005 Green
Description	Baby bonnet XS - Size 1 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 1 - 22-25 (cm) Head Circumference



Reference Code	AN0006 Yellow
Description	Baby bonnet S - Size 2 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 2 - 25-28 (cm) Head Circumference

Reference Code	AN0007 Red
Description	Baby bonnet M - Size 3 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 3 - 28-31 (cm) Head Circumference

> Patient circuit spare parts



Reference Code	AN0008 Blue
Description	Baby bonnet L - Size 4 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 4 - 31-34.5 (cm) Head Circumference



Reference Code	AN0009 Black
Description	Baby bonnet XL - Size 5 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 5 - 34.5-38 (cm) Head Circumference



Reference Code	016-012-000
Description	Silicone nasal prong
Size	size 0



Reference Code	016-012-001
Description	Silicone nasal prong
Size	size 1



Reference Code	016-012-002
Description	Silicone nasal prong
Size	size 2



Reference Code	016-012-003
Description	Silicone nasal prong
Size	size 3



Reference Code	016-012-004
Description	Silicone nasal prong
Size	size 4




Reference Code	016-012-005
Description	Silicone nasal prong
Size	size 5




Reference Code	016-012-006
Description	Silicone nasal prong
Size	size 6



Reference Code	016-012-007
Description	Silicone nasal prong
Size	size 7



Reference Code	016-001-000
Description	Nasal Cannula kit - Small
Includes	Nasal tubing, 1 nasal prongs (size 0), 1 baby bonnet, Velcro fixing strap and chin strap



Reference Code	016-001-001
Description	Nasal Cannula kit - Medium
Includes	Nasal tubing, 1 nasal prongs (size 1), 1 baby bonnet, Velcro fixing strap and chin strap

> Patient circuit spare parts



Reference Code	N4900 White
Description	RAM Cannula Micro Premie
Size	Size 0



Reference Code	N4901 Green
Description	RAM Cannula Premie
Size	Size 1



Reference Code	N4902 Blue
Description	RAM Cannula Newborn
Size	Size 2



Reference Code	N4903 Orange
Description	RAM Cannula Infant
Size	Size 3



Reference Code	N4904 Yellow
Description	RAM Cannula Small Paediatric
Size	Size 4



Reference Code	N4905 Teal
Description	RAM Cannula Medium Paediatric
Size	Size 5



Reference Code	N4906 Purple
Description	RAM Cannula Large Paediatric
Size	Size 6



Reference Code	1902000
Description	Y-piece for RAM Cannula
Quantity	1 / per patient per set-up



Reference Code	1666000
Description	Straight connector 15F - 10M Only required when using alongside autoclavable circuits (002-002-3150E)
Size	2 / per patient per set-up

> Water Trap (external) spare parts



Reference Code	1702000
Description	Elastomeric connector 15F - 6-9mm spigot
Quantity	1 / per water trap



Reference Code	1962000
Description	Straight connector 22F - 15M
Quantity	1 / per water trap



Reference Code	1997000
Description	Fixed elbow 15M-22M/15F
Quantity	1 / per water trap



Reference Code	1911000
Description	Water trap Self-sealing, medium volume
Quantity	1 / per water trap

Oxygen Concentrator spare parts <



Reference Code	AR20-F01B-B-CPAPMOD
Description	Regulator
Quantity	1



Reference Code	2510A2A14BVBN
Description	Acrylic Flow Meter - 0.4-5 LPM
Quantity	1



Reference Code	525DD-639
Description	Door intake filter package
Quantity	1



Reference Code	525DD-632
Description	Control circuit board starting with serial number B
Quantity	1



Reference Code	525D-619
Description	Sieve beds
Quantity	1



Reference Code	T25-01
Description	T25 screwdriver
Quantity	1

> Oxygen Concentrator spare parts




Reference Code	398 9091-MOD
Description	Water trap (Integral)
Quantity	1



Reference Code	812-106
Description	Straight connector - 6mm
Quantity	1



Reference Code	303DZ-605
Description	Cabinet Air Filter
Quantity	1



Reference Code	721-4007-A
Description	Black length tubing
Size	6mm x 1m



Reference Code	686-2680-A
Description	Clear length tubing
Quantity	6mm x 1m




Reference Code	748-0709
Description	Flow Controller
Size	4mm Inlet x 4mm Outlet

Oxygen Concentrator spare parts <



Reference Code	771-5370
Description	Push in straight fitting
Size	M5 x 6mm



Reference Code	91300000
Description	Voltage protector
Quantity	1



Reference Code	515A-705
Description	Exhaust muffler
Quantity	1



Reference Code	783-9574
Description	UK plug (Connector Converter)
Size	2M Length



Reference Code	42-0562
Description	EU (Power Lead)
Size	2.5M Length



Reference Code	42-0545
Description	UK (Power Lead)
Size	2M Length

> Oxygen Concentrator spare parts



Reference Code	505DZ-508
Description	Power switch
Quantity	1



Reference Code	MC44D-605
Description	Intake Bacteria Filter
Quantity	1



Reference Code	PV5LD-651
Description	Final Bacteria Filter
Quantity	1



Reference Code	515KZ-615
Description	Circuit Breaker, 3A
Quantity	1



Reference Code	501DZ-603
Description	Castor, non-locking
Quantity	Pack of 1



Reference Code	525DS-603
Description	Castor, Locking
Size	Pack of 1

Oxygen Concentrator spare parts <




Reference Code	525k-616
Description	Capacitor, 220V AC
Size	1



Reference Code	525K-625
Description	European compressor
Quantity	1



Reference Code	525DD-626
Description	Compressor Filter
Quantity	1



Reference Code	505IZ-609
Description	Motor Mounts
Quantity	Pack of 4



Reference Code	525K-643
Description	Compressor Rebuild Set
Quantity	1



Reference Code	515UK-634
Description	Cooling Fan
Quantity	1

> Oxygen Concentrator spare parts



Reference Code	525DD-614
Description	Intake Canister
Quantity	1



Reference Code	PVO2D-607
Description	Sieve Beds Check Valves
Quantity	1



Reference Code	515ADZ-702
Description	4-Way Rotary Valve
Quantity	1



Reference Code	525D-610
Description	Accumulator Tank
Quantity	1



Reference Code	MC29D-612
Description	Pressure Regulator
Quantity	1

BABY CPAP 20
SPARE PARTS LIST

> Baby CPAP 20 spare parts

CATEGORY	PART NAME	PAGE	CATEGORY	PART NAME	PAGE
Autoclavable circuits	Silicone autoclavable tubing	3	Water Trap (external) circuits	Elastomeric connector	10
	Straight connector	3		Straight connector	10
	Elastomeric connector	3		Fixed elbow - 15M-22M/15F	10
NON - Autoclavable circuits	Flextube	4		Water trap	10
	Elastomeric connector	4	Oxygen Concentrator circuits	Final bacteria filter	11
	Connector	4		Intake bacteria filter	11
	CPAP pipe fitting white acetal	4		Regulator	11
Patient circuit	Straight connector	4		Acrylic Flow Meter	11
	Humidifier bottle	5		Door intake filter package	11
	Bubble PAP valve bottle	5		Compressor service kit	11
	Baby bonnet XS - size 1	5		Control circuit board	12
	Baby bonnet S - size 2	5		Sieve beds	12
	Baby bonnet M - size 3	5		Capacitor	12
	Baby bonnet L - size 4	5		Exhaust muffler	12
	Baby bonnet XL - size 5	6		T25 screwdriver	12
	Silicone nasal prong - size 0	6		Water trap (Integral)	12
	Silicone nasal prong - size 1	6		Straight connector - 6mm	13
	Silicone nasal prong - size 2	6		Cabinet Air Filter	13
	Silicone nasal prong - size 3	6		Black length tubing	13
	Silicone nasal prong - size 4	6		Clear length tubing	13
	Silicone nasal prong - size 5	7		Voltage protector	13
	Silicone nasal prong - size 6	7		UK plug (Connector Converter)	13
	Silicone nasal prong - size 7	7		EU (Power Lead)	14
	Nasal Cannula kit - Small	7		UK (Power Lead)	14
	Nasal Cannula kit - Medium	7		Power switch	14
	RAM Cannula Micro Preemie - size 0	8		Compressor intake hose	14
	RAM Cannula Preemie - size 1	8		Compressor (230V) AC	14
	RAM Cannula Newborn - size 2	8		Cooling fan	14
	RAM Cannula Infant - size 3	8		Circuit breaker (5A)	15
	RAM Cannula Small Paediatric - size 4	8		Flow Controller	15
	RAM Cannula Medium Paediatric - size 5	8		Push in straight fitting	15
	RAM Cannula Large Paediatric	9		Castor, non-locking	15
	Y-piece for RAM Cannula	9		Castor, Locking	15
	Straight connector 15F - 10M	9		4-Way Rotary Valve	15
				1025D-610	16
				Pressure Regulator, 20PSI	16
				Motor Mounts	16

Contact information

Diamedica (UK) Ltd

Unit 1 Kingdom Fields,
Bratton Fleming,
Devon, UK
EX31 4EN

Tel: +44 (0)1598 710066

Email: info@diamedica.co.uk

Web: www.diamedica.co.uk

Technical Support:

Email: support@diamedica.co.uk

WhatsApp: +44 (0)7716 503156

Option one – Autoclavable Patient Limbs <



Reference Code	002-002-3150E
Description	Silicone autoclavable tubing - 10F cuff, length 150cm, reusable
Quantity	2 / per patient per set-up

Reference Code	1661000
Description	Straight connector - 15M/11.5 - 10M
Quantity	2 / per patient per set-up



Reference Code	1702000
Description	Elastomeric connector - 15F - 6-9mm spigot
Quantity	1 / per patient per set-up

> Option two – Single-use Patient Limbs



Reference Code	1576000
Description	Flextube - 15mm
Quantity	2 / per patient per set-up



Reference Code	1702000
Description	Elastomeric connector 15F - 6-9mm spigot
Quantity	1 / per patient per set-up



Reference Code	1945000
Description	Connector - 15m - 15m
Quantity	1 / per patient per set-up



Reference Code	PF-CPAP
Description	CPAP pipe fitting white acetal
Quantity	2 / per patient per set-up



Reference Code	1962000
Description	Straight connector - 22F - 15M
Quantity	1 / per patient per set-up

Patient circuit spare parts <



Reference Code	02-140
Description	7600-0 Humidifier bottle
Quantity	1



Reference Code	821135
Description	Reusable humidifier bottle
Quantity	1



Reference Code	BC01/AN0018
Description	Bubble PAP valve bottle
Quantity	1



Reference Code	AN0005 Green
Description	Baby bonnet XS - Size 1 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 1 - 22-25 (cm) Head Circumference



Reference Code	AN0006 Yellow
Description	Baby bonnet S - Size 2 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 2 - 25-28 (cm) Head Circumference



Reference Code	AN0007 Red
Description	Baby bonnet M - Size 3 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 3 - 28-31 (cm) Head Circumference

> Patient circuit spare parts



Reference Code	AN0008 Blue
Description	Baby bonnet L - Size 4 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 4 - 31-34.5 (cm) Head Circumference



Reference Code	AN0009 Black
Description	Baby bonnet XL - Size 5 Each Bonnet kit includes a bonnet, two adhesive hook style fasteners and ribbon.
Size	Size 5 - 34.5-38 (cm) Head Circumference



Reference Code	016-012-000
Description	Silicone nasal prong
Size	size 0



Reference Code	016-012-001
Description	Silicone nasal prong
Size	size 1



Reference Code	016-012-002
Description	Silicone nasal prong
Size	size 2



Reference Code	016-012-003
Description	Silicone nasal prong
Size	size 3

Patient circuit spare parts <



Reference Code	016-012-004
Description	Silicone nasal prong
Size	size 4




Reference Code	016-012-005
Description	Silicone nasal prong
Size	size 5




Reference Code	016-012-006
Description	Silicone nasal prong
Size	size 6



Reference Code	016-012-007
Description	Silicone nasal prong
Size	size 7



Reference Code	016-001-000
Description	Nasal Cannula kit - Small
Includes	Nasal tubing, 1 nasal prongs (size 0), 1 baby bonnet, Velcro fixing strap and chin strap



Reference Code	016-001-001
Description	Nasal Cannula kit - Medium
Includes	Nasal tubing, 1 nasal prongs (size 1), 1 baby bonnet, Velcro fixing strap and chin strap

> Patient circuit spare parts



Reference Code	N4900 White
Description	RAM Cannula Micro Preemie
Size	Size 0



Reference Code	N4901 Green
Description	RAM Cannula Preemie
Size	Size 1



Reference Code	N4902 Blue
Description	RAM Cannula Newborn
Size	Size 2



Reference Code	N4903 Orange
Description	RAM Cannula Infant
Size	Size 3



Reference Code	N4904 Yellow
Description	RAM Cannula Small Paediatric
Size	Size 4



Reference Code	N4905 Teal
Description	RAM Cannula Medium Paediatric
Size	Size 5

Patient circuit spare parts <



Reference Code	N4906 Purple
Description	RAM Cannula Large Paediatric
Size	Size 6



Reference Code	1902000
Description	Y-piece for RAM Cannula
Quantity	1 / per patient per set-up



Reference Code	1666000
Description	Straight connector 15F - 10M Only required when using alongside autoclavable circuits (002-002-3150E)
Size	2 / per patient per set-up

> Water Trap (external) spare parts



Reference Code	1702000
Description	Elastomeric connector 15F - 6-9mm spigot
Quantity	1 / per water trap



Reference Code	1962000
Description	Straight connector 22F - 15M
Quantity	1 / per water trap



Reference Code	1997000
Description	Fixed elbow 15M-22M/15F
Quantity	1 / per water trap



Reference Code	1911000
Description	Water trap Self-sealing, medium volume
Quantity	1 / per water trap

Oxygen Concentrator spare parts <



Reference Code	PV5LD-651
Description	Final bacteria filter
Quantity	1



Reference Code	1025D-605
Description	Intake bacteria filter
Quantity	1



Reference Code	AR20-F01B-B-CPAPMOD
Description	Regulator
Quantity	1



Reference Code	2510A2A15BVBN
Description	Acrylic Flow Meter 1-10 LPM
Quantity	1



Reference Code	525DD-639
Description	Door intake filter package
Quantity	1



Reference Code	1025K-643
Description	Compressor service kit
Quantity	1

> Oxygen Concentrator spare parts



Reference Code	1025D-622
Description	Control circuit board
Quantity	1



Reference Code	1025D-619
Description	Sieve beds
Quantity	1



Reference Code	1025UK-616
Description	Capacitor
Quantity	1



Reference Code	515A-705
Description	Exhaust muffler
Quantity	1



Reference Code	T25-01
Description	T25 screwdriver
Quantity	1



Reference Code	398 9091-MOD
Description	Water trap (Integral)
Quantity	1


Oxygen Concentrator spare parts <



Reference Code	812-106
Description	Straight connector - 6mm
Quantity	1



Reference Code	303DZ-605
Description	Cabinet Air Filter
Quantity	1



Reference Code	721-4007-A
Description	Black length tubing
Size	6mm x 1m



Reference Code	686-2680-A
Description	Clear length tubing
Quantity	6mm x 1m



Reference Code	91300000
Description	Voltage protector
Quantity	1



Reference Code	783-9574
Description	UK plug (Connector Converter)
Size	2M Length

> Oxygen Concentrator spare parts



Reference Code	42-0562
Description	EU (Power Lead)
Size	2.5M Length




Reference Code	42-0545
Description	UK (Power Lead)
Size	2M Length



Reference Code	1025D-508
Description	Power switch
Quantity	1



Reference Code	1025D-657
Description	Compressor intake hose
Quantity	1



Reference Code	1025K-625
Description	Compressor (230V) AC
Quantity	1



Reference Code	1025D-634
Description	Cooling fan
Quantity	1

Oxygen Concentrator spare parts <



Reference Code	1025K-613
Description	Circuit breaker (5A)
Quantity	1



Reference Code	748-0709
Description	Flow Controller
Size	4mm Inlet x 4mm Outlet



Reference Code	125-9596
Description	Push in straight fitting
Size	M5 x 6mm



Reference Code	501DZ-603
Description	Castor, non-locking
Quantity	Pack of 1



Reference Code	525DS-603
Description	Castor, Locking
Size	Pack of 1



Reference Code	515ADZ-702
Description	4-Way Rotary Valve
Quantity	1

> Oxygen Concentrator spare parts



Reference Code	1025D-610
Description	Accumulator Tank
Quantity	1



Reference Code	MC29D-612
Description	Pressure Regulator, 20PSI
Quantity	1



Reference Code	1025D-609
Description	Motor Mounts
Quantity	Pack of 4



Reference Code	PVO2D-607
Description	Beds Check Valves
Quantity	1